
BIOPOLYMER
CONNECT





We've obtained ISO 13485 Certification.

Quality has always been our main priority and our production is focused on your need.

We make an effort to provide products for our customers safely and reliably.
Also we are continuously working to provide superior service and innovative value-added medical solutions.

CERTIFICATION




AGECARE KOREA



*BIO POLYMER
CONNECT*

PRODUCT




GROA

Dermal Filler



VITA

PCL Threads



DIVE

Skin Booster




COMPLE

PLLA Powder

DOUBLE

An Innovative System for Skin Lifting
For Skin Firmness and Tension



BASIC

BASIC

04










COG MESH

COG

04

BASIC






High lifting effect can be expected due to COG's strong tensile strength and tissue fixation. It is effective for forehead and neck lifting. Depending on treatment area embedding can be used to maximize the effect.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT	-	2	320	○	-	-
	-	2	480	○	-	-
						
	-	2	320	-	○	○
	-	2	480	-	○	○
						
	-	2	320	○	○	-
	-	2	480	○	○	-

BLUNT 

COG MESH


High bearing capacity of Cog Mesh helps to achieve effective double chin lifting result. The procedure is simple and offers a long lasting effect.

NEEDLE		THREAD		MESH	PDO	PCL	PLLA	
Option	Gauge / mm	USP	mm	type	mm			
BLUNT	-	2	320	MOLD	220		-	-
	-	2	480	SLICE	380		-	-

BLUNT 

SINGLE

Stimulationg the Process of Neocollagenesis
Design Your Skin A Young Look



MESH

VOLUME MESH	06
CONE SPRING	06

COG

MOLD MESH	07
12D COG	07
VITA PCL MOLD	07
RIBBON MOLD	07
NOSE	08
COG SCREW	08












BASIC

MONO	09
MULTI MONO	09
MULTI PLLA 360	09
TWISTED SINGLE	10
TWISTED DOUBLE	10

MESH

VOLUME MESH














Multiple threads are braided for cells to move into the holes on the surface of maximizing the effect of tissue formation. Effective for removing deep-wrinkles.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	19G 38mm	6/0	40			-
BLUNT(A)	19G 50mm	6/0	60			-
BLUNT(L)	19G 60mm	6/0	80			-
BLUNT(L)	19G 100mm	6/0	160			-



CONE SPRING

Effective for removing wrinkles around nasolabial folds, marionette line.








NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	21G 60mm	2/0	80			-
	23G 60mm	3/0	160			-
	25G 50mm	5/0	40			-
SHARP	23G 60mm	3/0	60			
	21G 60mm	2/0	80			



COG

MOLD MESH






It has a strong lifting effect thanks to its traction and tensile strength due to its dual cutting.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	18G 60mm	5/0	80			-
BLUNT(L)	18G 100mm	5/0	160			-



12D COG

The thread has been processed similar to L-shape. The thread itself is thick so it is easy to anchor.







NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	19G 100mm	0	200			-



BEST PRODUCT

VITA PCL MOLD





Dual effect with Vitamin C & Polycaprolactone Copolymerization

NEEDLE		THREAD		PDO	V-PCL	TRI POLYMER
Option	Gauge / mm	USP	mm			
BLUNT(W)	18G 100mm	2	170			
BLUNT(L)	18G 100mm	2	170			



RIBBON MOLD




Double sided cog catches and pulls the tissue well to give a natural lifting effect.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(W)	19G 100mm	6/0	160	-	-	






NOSE

It is optimized for cosmetology.
In 19G 38mm it was used for
Columella, in 19G 50mm for
dorsum.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(H-L)	19G 38mm	0	40	○	○	-

BLUNT (H-L) 






NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(H-L)	19G 50mm	2	70	○	○	-

BLUNT (H-L) 



COG SCREW

It reaches deeper level of the skin
and covers larger area.
Therefore, tissue stimulation effect
is higher.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	23G 60mm	3/0	195	-	-	○

BASIC 



BASIC

BEST PRODUCT

MONO

HA PCL - Dual effect with Sodium Hyaluronate & Polycaprolactone Copolymerization

NEEDLE		THREAD		PDO	HA PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	26G 60mm	5/0	90	○	○	○
BASIC	27G 60mm	5/0	90	○	○	○
BASIC	29G 38mm	6/0	50	○	○	○
BLUNT(F)	30G 25mm	6/0	30	○	○	-



MULTI MONO

Six mono thread can be in inserted at once, effective for neck wrinkles.


NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	21G 60mm	5/0	90	○	○	○



MULTI PLLA 360

360 mono thread can be in inserted at once, effective for neck wrinkles.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	23G 50mm	9/0	60	-	-	○
BASIC	27G 38mm	11/0	50	-	-	○



TWISTED SINGLE

Excellent for tightening effect.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	26G 38mm	5/0	50	○	○	○
BASIC	26G 60mm	5/0	90	○	○	○
BASIC	29G 38mm	6/0	50	○	○	○
BASIC	29G 60mm	6/0	90	○	○	○
BASIC	30G 25mm	6/0	30	○	○	○

BASIC



TWISTED DOUBLE

Commonly used for marionette line, short procedure time.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	26G 38mm	7/0	50	○	○	○
BASIC	26G 60mm	7/0	90	○	○	○

BASIC





#816 31, Bongmyeong-ro, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, Republic of Korea
TEL +82 10 5481 0512 **FAX** +82 0504 162 2362 **E-MAIL** annaleeack@gmail.com
