
LIQUID
THREAD
REAL

BIOMEDICAL
MATERIALS

DEVICES
ACROSS
TECHNOLOGY




**AGECARE
KOREA**



DOUBLE

An Innovative System for Skin Lifting
For Skin Firmness and Tension



BASIC

BASIC04

COG MESH

COG04

BASIC

High lifting effect can be expected due to COG's strong tensile strength and tissue fixation. It is effective for forehead and neck lifting. Depending on treatment area embedding can be used to maximize the effect.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT	-	0	290		-	-
	-	0	320		-	-
	-	2/0	440		-	-
	-	2/0	480		-	-

BLUNT

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT	-	0	290	-		
	-	0	320	-		
	-	2/0	440	-		
	-	2/0	480	-		

BLUNT

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT	-	0	290			-
	-	0	320			-
	-	2/0	440			-
	-	2/0	480			-

BLUNT

COG


High bearing capacity of Flat Mesh helps to achieve effective double chin lifting result. The procedure is simple and offers a long lasting effect.

NEEDLE		THREAD		MESH	PDO	PCL	PLLA	
Option	Gauge / mm	USP	mm	type	mm			
BLUNT	-	0	145×2	MOLD	100		-	-
	-	1	145×2	SLICE	100		-	-
	-	2	145×2	SLIEC	100		-	-

BLUNT

SINGLE

Stimulationg the Process of Neocollagenesis
Design Your Skin A Young Look



MESH

VOLUME MESH	06	DUAL COG	07	MONO	09
CONE MESH	06	12D COG	07	MULTI MONO	09
		MOLD COG	07	MULTI PLLA 360	09
		RIBBON COG	07	TWISTED SINGLE	10
		NOSE	08	TWISTED TWIN	10
		COG SCREW	08	TWISTED DOUBLE	10
				TRIPLE DOUBLE	10

MESH

VOLUME MESH

Multiple threads are braided for cells to move into the holes on the surface of maximizing the effect of tissue formation. Effective for removing deep-wrinkles.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	18G 60mm	5/0	80			-
BLUNT(L)	18G 100mm	5/0	160			-
BLUNT(L)	19G 38mm	6/0	40			-
BLUNT(A)	19G 50mm	6/0	60			-
BLUNT(L)	19G 60mm	6/0	80			-
BLUNT(L)	19G 100mm	6/0	160			-



CONE MESH

Effective for removing wrinkles around nasolabial folds, marionette line.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	18G 60mm	5/0	80			-
BLUNT(L)	18G 100mm	5/0	160			-
BLUNT(L)	19G 38mm	6/0	40			-
BLUNT(A)	19G 50mm	6/0	60			
BLUNT(L)	19G 60mm	6/0	80			
BLUNT(L)	19G 100mm	6/0	160			-



COG

DUAL COG

It has a strong lifting effect thanks to its traction and tesile strength due to its dual cutting.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	18G 60mm	5/0	80			-
BLUNT(L)	18G 100mm	5/0	160			-
BLUNT(L)	19G 38mm	6/0	40			-
BLUNT(A)	19G 50mm	6/0	60			-
BLUNT(L)	19G 60mm	6/0	80			-
BLUNT(L)	19G 100mm	6/0	160			-



12D COG

The thread has been processed similar to L-shape. The thread itself is thick so it is easy to anchor.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	19G 100mm	0	200			-



MOLD COG

It is effective for deep and dark wrinkles. Moreover it has strong lifting effect.

NEEDLE		THREAD		PDO	PCL	TRI POLYMER
Option	Gauge / mm	USP	mm			
BLUNT(W)	18G 100mm	2	170			
BLUNT(L)	18G 100mm	2	170			



RIBBON COG

Double sided cog catches and pulls the tissue well to give a natural lifting effect.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(W)	19G 100mm	6/0	160	-	-	



NOSE

It is optimized for cosmetology. In 19G 38mm it was used for Columella, in 19G 50mm for dorsum.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(H-L)	19G 38mm	0	40			-



NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(H-L)	19G 50mm	2	70			-



COG SCREW

It reaches deeper level of the skin and covers larger area. Therefore, tissue stimulation effect is higher.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	25G 60mm	5/0	195	-	-	



BASIC

MONO

It gives natural volume and used as a pain treatment.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	26G 38mm	5/0	50	○	○	○
BASIC	26G 60mm	5/0	90	○	○	○
BASIC	27G 38mm	5/0	50	○	○	○
BASIC	27G 60mm	5/0	90	○	○	○
BASIC	29G 38mm	6/0	50	○	○	○
BASIC	29G 50mm	6/0	80	○	○	○
BASIC	29G 60mm	6/0	90	○	○	○
BLUNT(F)	30G 25mm	6/0	30	○	○	-



MULTI MONO

Six mono thread can be in inserted at once, effective for neck wrinkles.


NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	21G 60mm	5/0	90	○	○	○



MULTI PLLA 360

Six mono thread can be in inserted at once, effective for neck wrinkles.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BLUNT(L)	23G 50mm	9/0	60	-	-	○
BASIC	27G 38mm	11/0	50	-	-	○



TWISTED SINGLE

Excellent for tightening effect.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	26G 38mm	5/0	50	○	○	○
BASIC	26G 60mm	5/0	90	○	○	○
BASIC	29G 38mm	6/0	50	○	○	○
BASIC	29G 60mm	6/0	90	○	○	○
BASIC	30G 25mm	6/0	30	○	○	○



TWISTED TWIN

Two threads are twisted form for maximizing tissue stimulation, used for face and neck wrinkles.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	26G 38mm	7/0	50	○	○	○
BASIC	26G 60mm	7/0	80	○	○	○
BASIC	36G 90mm	7/0	130	○	○	○



TWISTED DOUBLE

Commonly used for marionette line, short procedure time.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	26G 38mm	7/0	50	○	○	○
BASIC	26G 60mm	7/0	90	○	○	○
BASIC	36G 90mm	7/0	130	○	○	○



TRIPLE

Effective for fine lines.

NEEDLE		THREAD		PDO	PCL	PLLA
Option	Gauge / mm	USP	mm			
BASIC	25G 38mm	7/0	50	○	○	○
BASIC	25G 60mm	7/0	90	○	○	○





#816 31, Bongmyeong-ro, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, Republic of Korea
TEL +82 10 5481 0512 **FAX** +82 0504 162 2362 **E-MAIL** annaleeack@gmail.com
